
HD 3.5 Portable Digital Video Magnifier

User's Guide

Content

Welcome To Use.....	2
Product Features.....	2
Packaging & Accessories.....	3
Technical Specifications.....	4
Product Structure.....	5
Use Introduction.....	6
1、 Install battery.....	6
2、 Battery charging.....	6
3、 Power on / off.....	7
4、 Use digital video magnifier	7
Functions Introduction.....	8
1、 Color model switch.....	8
2、 Zoom rate adjustment.....	8
3、 Image freeze/unfreeze.....	9
4、 Reading line.....	9
5、 Flashlight.....	9
6、 Memory function.....	10
7、 Power saving.....	10
8、 Low battery tips	10
9、 Handwriting bracket.....	11
10、 Connecting to TV.....	11
Menu Introduction.....	12
Computer configuration requirements.....	13
Precautions.....	13
Problem solutions.....	14

Welcome To Use

Thank you for purchasing our company's portable digital video magnifier! Please read and keep this manual carefully. It will help you to use the magnifier better and achieve the best magnifying effect.

Product Features

This Full HD 3.5 inch portable digital video magnifier is our company grand launch with high performance and high cost-effective performance.

- 720P60 design, no lag, no smear, fluent HD.
- Professional image processing DSP to restore colors and details really;
- Intelligent image enhancement algorithm; Adapting to image scene intelligently and reproducing details after zooming in;
- 3.5 inch large anti-glare screen, 960X240 resolution, high definition and high lighting;
- Custom-made high-end lens, large filed angle and extended depth of field, more text and continuous reading;
- 2X-32X infinite amplification, clear and smooth;
- Ultrathin style and portable; Zero distance design, use comfortable;
- Concave two-color button, high contrast, clear touch;
- Ergonomic bracket design, comfort to reading and fluent handwriting;
- 26 color modes to meet different people's need;
- Color mode custom-made function, any combination of enhancing mode, personalized custom;
- Picture freeze, convenient to take a close look;
- Reading line function and auxiliary continuous reading text;
- Large icon user interface, clear and concise;
- Adjustable backlight, comfortable, eye protected; White LED light brightness can be adjusted or turn off, Suitable for reading mobile phone, tablets and other reading materials;
- Memory function can automatically save the user all the last settings and remember the user's habits;

-
- Support AV and connecting to HD TV to provide greater and clearer visual effects;
 - Support external headphones to avoid the prompt tone affecting the surrounding environment
 - Flashlight can supply temporary lighting;
 - Large capacity removable battery, the ultimate power-saving design, 4 hours runtime;
 - Intelligent power-saving design; Automatic shut-down after the screen is static and no operation for 3 minutes ; Environmental protection and energy conservation;

New algorithm, new design and new configuration will bring you smooth, high-definition, portable new experience!






Packaging & Accessories

- | | | | |
|---------------------------|-----|-----------------------|-----|
| ● Digital video magnifier | 1PC | ● User manual | 1PC |
| ● Video output cable | 1PC | ● Sling | 1PC |
| ● USB charging cable | 1PC | ● Protective case | 1PC |
| ● Power adapter | 1PC | ● Lens cleaning cloth | 1PC |
| ● Lithium battery | 1PC | | |

Technical Specifications

- Screen size: 3.5 HD colorful LCD screen(960X240)
- Pixel size: 1 Mega pixel (720P60)
- Zoom rate: 2X ~ 32X Infinite amplifying
- Color mode: 26 kinds
- Mode custom-made: Support
- White LED adjustment: Support
- Brightness adjustment: Support
- Memory function: Support
- Image freeze: Support
- Reading line: Support
- Flashlight: Support
- TV output: Support AV
- Headphone output: Support 2.5mm headphone
- Power saving function: Support
- Low battery tips: Support
- Voice prompt: Support
- Handwriting bracket: Support
- Working time: last over 4 hours
- Battery capacity: 1500mAH high-capacity chargeable lithium battery
- Dimension: 115mm(length)x 75mm(width)x 18mm(height)
- Weight: 150g(including the battery)


Product Structure

- | | | |
|---|----------------------------|-------------------------|
| 1. Power button  | 7. Charging indicator(red) | 13. Lens gobe |
| 2. Freeze button  | 8. LCD gobe | 14. White LED |
| 3. Zoom out button  | 9. Lens focus | 15. Handwriting bracket |
| 4. Zoom in button  | 10. AV interface | 16. Battery cover |
| 5. Mode button M | 11. USB interface | 17. Speaker |
| 6. Power indicator(green) | 12. Wrist strap shaft | 18. Battery |


Use Introduction

1. Install battery

- Invert video magnifier;
- Hold the battery cover and push the cover towards edge to open;
- Install the battery according to the positive and negative tips;
- Cover the battery box, sliding in place when hearing the "click" sound.

2. Battery charging

- Connecting video magnifier to adapter by USB cable;
- Connecting adapter to AC power;
- Charging light(red light) is on when charging;
- Charging light(red light) will be off when complete charge;
- At the first time, please charge the battery to full
- When low battery tips, please charge the battery in time

Attention:

In the charging process, please keep the charger and the video magnifier have good heat dissipation and ventilation, so as to avoid danger caused by overheat.

3. Power on/off

Long press power button to switch on

Long press power button again to switch off

- Keep pressing **power button**  for 3 seconds while the power is off, the video magnifier will turn on, and the green indicator lights up, and video magnifier enters the vision aid mode.
- Aim the lens focus at the reading target;
- Optimal reading effect can be attained through adjusting magnifying power, enhanced mode, locking picture and the like operations;
- Keep pressing **power button**  for 3 seconds while the power is on, the video magnifier will turn off;

Attention:

For best results, remove the protective film from the LCD and lens goggles before using.

4. Use video magnifier

Color mode switch: Press **Mode button**  to switch color mode;

Zoom in and out : Press **Zoom in button**  to enlarge image and press **Zoom out button**  to shrink image.

Freeze and unfreeze image: Press **freeze button**  to freeze image and repress **freeze button**  to unfreeze image.

System parameter setting: Long press **Mode button**  to enter menu setting and long press **freeze button**  to enter information menu.

For more details, please refer to the chapter on function introduction and menu setting.

Functions Introduction

1. Color mode switch


Switch the LCD screen image enhancement mode by Mode button to meet different needs.

26 kinds of image enhancement modes:

- | | |
|--------------------------|--------------------|
| 1)Color mode | 14)White on green |
| 2)Color enhancement mode | 15)Orange on white |
| 3)Negative mode | 16)White on orange |
| 4)Monochrome mode | 17)Red on black |
| 5)Black on white | 18)Black on red |
| 6)White on black | 19)Orange on black |
| 7)Blue on white | 20)Black on orange |
| 8)White on blue | 21)Yellow on green |
| 9)Yellow on black | 22)Green on yellow |
| 10)Black on yellow | 23)Red on white |
| 11)Yellow on blue | 24)White on red |
| 12)Blue on yellow | 25)Yellow on white |
| 13)Green on white | 26)White on yellow |

The desired color mode can be customized by menu setting -> color mode, details see Color mode setting section.

2.Zoom rate adjustment

2X Amplifying

16X Amplifying

32X Amplifying

By zoom in or zoom out to adjust zoom rate to achieve the best magnifying effect
This magnifier can zoom in or zoom out infinitely between 2X ~ 32X.

3. Image freeze and unfreeze

Freeze image after viewfinder


View image after freezing

Freeze button  can freeze image and can observe the image at close distance, which is not easy to watch before freezing;


Repress **Freeze button**  to unfreeze image and recover real-time view.

4. Reading line


Reading lines can help text positioning to help users read quickly and continuously.

Horizontal reading line


Vertical reading line

Turn off the reading line

Long press **Freeze button**  to switch horizontal reading line, vertical reading line and close the reading line.

Moving reading line : When the reading line is flashing, the reading line is moved by pressing **Zoom in**  and **Zoom out**  buttons.

5. Flashlight

Flashlight plays the role of temporary lighting for video magnifier.

Turn on flashlight: under off status, keep pressing **freeze button**  and then press power button, flashlight is on.

Mode switch: By **freeze button**  to switch bright, slow flash, flash, gradient and close.

Brightness adjustment: **Zoom in +** or **zoom out -** can brighten or dim flashlight brightness (except for the gradient mode)

Turn off flashlight: Long press power button and flashlight will be off.

6.Memory function

After turning off, memory function can save the user all the last settings automatically and remember the user's habits;

7.Power saving


The video magnifier has intelligent power saving function. When no key operation and the screen rests for 1 minute, it will enter the sleep state automatically and the screen shows the battery power and the current date and time.

Automatically sleep

Under sleep status, pressing any key to wake up magnifier. No operation after 2 minutes, it will shut down automatically to save power.

8.Low battery tips

Low battery tips

Low battery tips function: Screen will display the low power icon at low battery, and the power light flashes.

When low battery tips, please charge the battery in time.

9. Handwriting bracket

This magnifier supports handwriting function. Open the bracket to write.


Writing on bracket

Rotate the magnifier 180 degrees for left hand writing.

10. Connecting to TV

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Connect magnifier to TV by AV video cable to get a bigger, clearer display.

When the video magnifier enters into video output mode, the LCD screen will turn off automatically.

Menu Introduction

Long press **Mode button M** to open Menu and set the parameters;

Under menu open mode, the buttons are defined as follows:

Pressing button	Operation
Mode button M Short press	Next menu
Zoom in +	Previous menu
Freeze button 	Option setting
Zoom out -	Next option
Mode button M Long press	Exit menu

Under Preview Mode, long press **Mode button M** to enter Menu and set the parameters of the video magnifier:

LCD brightness adjustment: can adjust between 1 ~ 5.

White LED brightness adjustment: can adjust between 0 ~ 5, 0 level to turn off white LED, mainly for reading mobile phones, panel and texts on other self-luminous instruments.

Beep volume adjustment: can adjust between 0 ~ 5, 0 level is beep off .

Date and time setting: can set current date and time;

Color mode setting: can choose the required using color mode, color mode cannot be turned off.

Battery power display: Show the current battery level as a percentage;

Version information display: Show the product model, firmware version number, and firmware compilation time.

Computer configuration requirements

This digital video magnifier requires the PC to meet the following requirements:

- Windows® Win7/Vista/XP operating system;
- Intel® Pentium 4, 2.8 GHz CPU or above;
- At least 1GB RAM;
- Standard USB 2.0 or above;
- 512MB or above video card;

Precautions:

In order to keep digital video magnifier works in the best condition, please read the following notes carefully:

- Please do not expose it directly to high temperatures or sunlight;
- Please keep it away from moisture, rain, liquids and chemicals;
- Please do not use it near to electronic equipments, and medical equipments, which are not under insufficient shield;
- Please ensure that using it at temperatures between 10 ° C and 40 ° C;
- Please ensure that it is stored at temperatures between -20 ° C and 65 ° C;
- Do not try to repair or open the aids to avoid damaging warranty;
- Please use the matching power supply and battery, otherwise danger or loss will occur;
- Unplug power and remove battery before cleaning magnifier. Use soft cloth to clean it, alcohol solvents or any other cleaners cannot be used

Warning:

Don't use this magnifier during charging.

Problem solutions

Problem	Solutions
Black screen	Battery is exhausted. Please recharge. Make sure the lens focus on the reading.
cannot power on	Battery is exhausted. Please recharge.
No sound	Confirm whether the volume at 1 to 5 in the menu setting
The image is too dark or too bright.	Make sure the lens focus on the reading.
Image is blurred or dirty	Wipe the LCD screen and lens goggles.
Reading text is too small.	To adjust zoom rate by pressing the zoom in button.
Connected to the TV, but the TV is not shown.	Make sure the TV is on. When using the AV cable, switch the TV input to AV input. Make sure the video cable is properly connected to the magnifier and TV.

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