

Record Weight Option (continued):

3. To record your weight, repeat steps 1-2. Then, press ▲ button (F) to select USER 1 or 2 (C). Stand on scale steadily. Your weight will display and then read out loud twice. Weight data will be saved in chronological order. Each user can save their weight up to 30 times (A).
- To view saved weight record, press MODE button (E). Then press ▲ button (F) to select USER 1 or 2 (C). Press MODE button (E) to display the latest weight record (B). The word "TIME" (A) will flash. Press ▲ button (F) or ▼/U button (G) to view 1- 30 recorded weight records (A). Note: weight records will not read out loud.
 - To delete saved weight record, press MODE button (E). Then press ▲ button (F) to select USER 1 or 2 (C). Then press MODE button (E) and release. Then press MODE button (E) and ▼/U button (G) at the same time for a few seconds. "del" will appear on display and weight record will reset to 0,0. Repeat to delete each USER's weight record.

Trouble Shooting:

- "LO" means low battery. Replace with new battery.
- "OL" means scale is overloaded.

Cleaning:

- Clean surface with a cloth. Never use solvents or abrasive cleaning agents.



For additional assistance or trouble-shooting help,
please email info@jobar.com

© 2018 Jobar International, Inc.
Carson, CA 90745
Made in China

**Extra Wide Talking Scale
Instructions**

IMPORTANT: PLEASE READ ALL THE INSTRUCTIONS CAREFULLY BEFORE USING THIS PRODUCT.

Battery Precautions:

- If battery contents are leaking, discontinue use and discard.
- Remove dead batteries or if not in use.
- Make sure to properly dispose dead batteries in accordance with state and local laws.

Battery Installation:

- Slide battery cover (H) to open and insert a 9V battery (not included). Replace battery cover (H). Note: battery compartment is design to fix battery snug.

Scale Functions:

Note: for accurate weight reading do not place scale on top of carpet floors. Make sure to place scale on flat surface only.

1. Step on scale platform or press MODE button (E) to turn ON.
2. Press ▼/U button (G) to select weight in lbs or kg. Then stand on scale steadily. Your weight will display and then read out loud twice.

Weight Specification:

- Max weight: 550 lbs (250 kg).

